Cat. No. BME100749



PRODUCT INFORMATION

Uniprot ID P59901

HZN-7734;HZN7734;MEDI-7734;MEDI7734;VIB-7734;VIB7734 **Common Name**

Conjugate Unconjugated ILT7;CD85g **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of Formulation & Reconstitution

reconstitution.

Host Species lgG1 IgG type Reactivity Human LILRA4 **Target**

Description Anti-LILRA4(daxdilimab biosimilar∏lgG1) mAb

Delivery 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Usage Research use only

