

## **PRODUCT INFORMATION**

**Uniprot ID** Q8NHJ6 TBC **Common Name** 

Conjugate Unconjugated

**Synonyms** B4;ILT3;LIR5;CD85K;ILT-3;LIR-5;HM18

**Applications** N/A

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution for specific instructions of reconstitution.

**Host Species** N/A IgG4 IgG type Reactivity Human **Target** LILRB4

**Description** Anti-LILRB4(rolistobart biosimilar IgG4) mAb

**Delivery** 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

