

PRODUCT INFORMATION

Common Name	PR-1498487
Conjugate	Unconjugated
Synonyms	LIB, LRC15
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Chimeric Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	LRRC15
Uniprot ID	Q8TF66
Description	Anti-LRRC15(samrotamab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

CONFIDENTIAL






Figure 1. Flow cytometry analysis with 1 μ g/mL Anti-LRRC15(samrotamab biosimilar) mAb (BME100111) on HEK293 cells transfected with Human LRRC15 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

DIMA BIO CONFIDENTIAL

