Cat. No. BME100623



PRODUCT INFORMATION

Uniprot ID P10636

JNJ-63733657;JNJ 63733657;JNJ63733657;JNJ-3657;JNJ3657;JNJ 3657 **Common Name**

Coniugate Unconjugated

TAU;FTD1;MSTD;PPND;DDPAC;MAPTL;MTBT1;MTBT2;tau-40;FTDP-17;PPP1R103;Tau-PHF6 Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type lgG1 Reactivity Human MAPT Target

Anti-MAPT(posdinemab biosimilar∏lgG1) mAb Description

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Background

Research use only

Usage



