

## **PRODUCT INFORMATION**

**Uniprot ID** P10636

ABBV-8E12;ABBV8E12;ABBV 8E12 **Common Name** 

Coniugate Unconjugated

TAU;FTD1;MSTD;PPND;DDPAC;MAPTL;MTBT1;MTBT2;tau-40;FTDP-17;PPP1R103;Tau-PHF6 Synonyms

**Applications** Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** IgG type lgG4 Reactivity Human MAPT Target

Anti-MAPT(tilavonemab biosimilar∏lgG4) mAb Description

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Background 

Usage Research use only

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

