

PRODUCT INFORMATION

Clone ID	OTI4E10 (formerly 4E10)
Target	mOrange, mOrange2
Synonyms	N/A
Host Species	Mouse
Description	Carrier-free (BSA/glycerol-free) Mouse monoclonal mOrange/mOrange2 antibody, clone OTI4E10
Delivery	3-4week
Uniprot ID	N/A
IgG type	IgG2b
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:2000
Purification	N/A
Formulation & Reconstitution	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

