

PRODUCT INFORMATION

Uniprot ID P23510

Common Name KY-1005;KY1005;SAR-445229

Conjugate Unconjugated

GP34;CD252;OX4OL;TXGP1;CD134L;OX-40L;TNLG2B;TNFSF4;Glycoprotein Gp34;CD134 ligand **Synonyms**

Applications Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG4 IgG type Reactivity Human OX40 Ligand Target

Description Anti-OX40 Ligand(amlitelimab biosimilar ☐ IgG4) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are Storage & Shipping

shipped at ambient temperature.

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Usage Research use only