

PRODUCT INFORMATION

Uniprot ID Q8NBP7

SHR-1209;SHR1209;SHR 1209 **Common Name**

Conjugate Unconjugated

Synonyms FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human **Target** PCSK9

Description Anti-PCSK9(recaticimab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

