

PRODUCT INFORMATION

Uniprot ID 015116

Common Name Libtayo; REGN-2810; REGN2810; SAR-439684; SAR439684

Conjugate Unconjugated

Synonyms PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific Formulation &

Reconstitution

instructions of reconstitution.

Host Species N/A IgG type laG4 Reactivity Human PD-1 **Target**

Description Anti-PD-1(cemiplimab biosimilar ☐ lgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped

at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only **Usage**

