

## **PRODUCT INFORMATION**

Uniprot ID Q15116

Common Name APL 501;APL-501;APL501;CBT

501;CBT501;CBT-501;GB 226;GB-226;GB226

Conjugate Unconjugated

**Synonyms** PD1;PD-1;CD279;SLEB2;hPD-1;hPD-l;hSLE1

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG4
Reactivity Human
Target PD-1

**Description** Anti-PD-1(geptanolimab biosimilar IgG4) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

**Storage & Shipping** intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in

Email: info@dimabio.com Website: www.dimabio.com

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

