Cat. No. BME100609



## **PRODUCT INFORMATION**

**Uniprot ID** Q9NZQ7

**Common Name** CX-072;CX072;CX 072

Conjugate Unconjugated

B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1 **Synonyms** 

**Applications** N/A Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** IgG type lgG4 Reactivity Human PD-L1 **Target** 

**Description** Anti-PD-L1(pacmilimab biosimilar IgG4) mAb

**Delivery** 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

