Cat. No. BME100679



PRODUCT INFORMATION

Uniprot ID Q9NZQ7

Common Name HS636;HS-636;HS 636

Conjugate Unconjugated

B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1 **Synonyms**

Applications N/A Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type lgG1 Reactivity Human PD-L1 **Target**

Description Anti-PD-L1(sudubrilimab biosimilar IgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

