Cat. No. BME100751



PRODUCT INFORMATION

Uniprot ID Q96B86

ABT-555;ABT555 **Common Name** Conjugate Unconjugated

Synonyms RGM Applications N/A

Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human **Target RGMA**

Description Anti-RGMA(elezanumab biosimilar□lgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

