

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | Q02161 |
| Common Name | r-anti-D |
| Conjugate | Unconjugated |
| Synonyms | RH;Rh4;RH30;RhII;RhPI;DIIIc;RHCED;RHDeI;RHPII;RhDCw;CD240D;HDFNRH;RHXIII;SLC42A5;RHDVA(TT);RhK562-II |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | RHD |
| Description | Anti-RHD(trinbelimab biosimilar) IgG1 mAb |
| Delivery | 4-6 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

