

PRODUCT INFORMATION

Uniprot ID Q6ZMC9

Common Name

BSI 060T;BSI-060T;BSI060T;PYX

106;PYX-106;PYX106

Conjugate Unconjugated

Synonyms CD33L3;HsT1361;SIGLEC-15

Applications N/A
Recommended N/A
Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target SIGLEC15

Description Anti-SIGLEC15(aldastotug biosimilar [lgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Storage & Shipping intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

