

## **PRODUCT INFORMATION**

**Uniprot ID** P78324

BMS-986351;BMS986351;CC-95251;CC95251 **Common Name** 

Conjugate Unconjugated

**Synonyms** BIT;MFR;P84;SIRP;MYD-1;SHPS1;CD172A;PTPNS1

**Applications** N/A

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** N/A lgG1 IgG type Reactivity Human **Target** SIRPα

**Description** Anti-SIRPα(anzurstobart biosimilar□lgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

