

## **PRODUCT INFORMATION**

**Uniprot ID** Q9BQB4

BPS-804;BPS804;BPS

804;MOR-05813;MOR05813;MOR **Common Name** 

05813;MOR-5813;MOR5813;MOR 5813

Conjugate Unconjugated

CDD;VBCH;DAND6;SOST1 **Synonyms** 

**Applications** N/A Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

**Host Species** N/A IgG type IgG2 Reactivity Human **SOST Target** 

**Description** Anti-SOST(setrusumab biosimilar IgG2) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background** 

therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

or animals.

**Usage** Research use only

