

PRODUCT INFORMATION

Clone ID	OTI2H8 (formerly 2H8)
Target	tGFP
Synonyms	N/A
Host Species	Mouse
Description	Mouse monoclonal turboGFP antibody, clone OTI2H8
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG2b
Clonality	Monoclonal
Reactivity	N/A
Applications	WB IF
Recommended Dilutions	WB: 1:1000-1:2000, IF: 1:100
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage & Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated




Figure 1. Immunoblot analysis of randomly picked 7 different tGFP tagged over-expression lysates with 2H8 anti-tGFP monoclonal antibody (DME900005) at 1:2000 dilution.






Figure 2. Immunoprecipitation (IP) of recombinant tGFP fusion protein from overexpression lysate made with HEK293T cells transfected by plasmid expressing tGFP-DDK. The immunoprecipitated tGFP fusion protein was analyzed by immunoblot with rabbit anti-DDK polyclonal antibody. (Left lane: Negative control, IP without adding DME900005)

DIMABIO CONFIDENTIAL

