

**PRODUCT INFORMATION**

<b>Clone ID</b>	OT12H8 (formerly 2H8)
<b>Target</b>	tGFP
<b>Synonyms</b>	N/A
<b>Host Species</b>	Mouse
<b>Description</b>	Mouse monoclonal turboGFP antibody, clone OT12H8
<b>Delivery</b>	1-2week
<b>Uniprot ID</b>	N/A
<b>IgG type</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	N/A
<b>Applications</b>	WB IF
<b>Recommended Dilutions</b>	WB: 1:1000-1:2000, IF: 1:100
<b>Purification</b>	N/A
<b>Formulation &amp; Reconstitution</b>	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Storage&amp;Shipping</b>	Store at -20°C.
<b>Background</b>	N/A
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

DIMABIO CONFIDENTIAL

