

## **PRODUCT INFORMATION**

Uniprot ID Q495A1

Common Name

BGB-A1217;BGB A1217;WCD-118;WCD118;WCD

118

**Conjugate** Unconjugated

Synonyms VSIG9;VSTM3;WUCAM

Applications N/A Recommended

Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target TIGIT

**Description** Anti-TIGIT(ociperlimab biosimilar □lgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

**Storage & Shipping** intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only



