


PRODUCT INFORMATION


Common Name	TSR-022
Synonyms	HAVCR2;TIM3;TIMD3;FLJ14428;KIM3
Conjugate	Unconjugated
Applications	ELISA; Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000; Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG4 - Kappa
Reactivity	Human
Target	TIM3
Uniprot ID	Q8TDQ0
Description	Anti-TIM3 (cobolimab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Usage	Research use only
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.



**Anti-TIM3 (cobolimab biosimilar) mAb ELISA**  
0.2 µg of Human TIM3, mFc-His Tagged protein per well



**Figure 1.** ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human TIM3, mFc-His tagged protein ([getskuurl sku="PME100030"]) can bind Anti-TIM3 Neutralizing antibody in a linear range of 0.64-80 ng/ml.



**Figure 2.** Flow cytometry analysis with 1µg/mL Anti-TIM3 (cobolimab biosimilar) mAb (BME100008) on HEK293 cells transfected with Human TIM3 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

