

PRODUCT INFORMATION

Uniprot ID Q8TDQ0

Common Name BGB-A425;BGB A425

Conjugate Unconjugated

Synonyms TIM3;CD366;KIM-3;SPTCL;TIMD3;Tim-3;TIMD-3;HAVcr-2

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % -8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

Host Species IgG type lgG1 Reactivity Human **Target** TIM3

Anti-TIM3(surzebiclimab biosimilar IgG1) mAb Description

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

