Cat. No. BME100724



PRODUCT INFORMATION

Uniprot ID Q8TDQ0

Common Name INCAGN2390;INCAGN02390;INCAGN-2390;INCAGN-02390

Coniugate Unconjugated

Synonyms TIM3;CD366;KIM-3;SPTCL;TIMD3;Tim-3;TIMD-3;HAVcr-2

Applications Recommended N/A Dilutions

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions

of reconstitution.

Host Species IgG type lgG1 Reactivity Human **Target** TIM3

Description Anti-TIM3(verzistobart biosimilar [IgG1) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping freezing and thawing). Lyophilized proteins are shipped

at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Research use only Usage

