PRODUCT INFORMATION

Uniprot ID 000206

Common Name EB-05, NI-0101 Conjugate Unconjugated

toll like receptor 4, toll-like receptor 4, hToll, CD284, TLR-4 **Synonyms**

Applications Flow Cyt

Recommended

Flow Cyt 1:100 Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species

IgG type Human IgG1(K214R,N325S,L328F) - Kappa

Reactivity Human TLR4 **Target**

Description Anti-TLR4(paridiprubart biosimilar) mAb

Delivery In Stock

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized antibodies are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only Usage








Figure 1. Flow cytometry analysis with 1 μ g/mL Anti-TLR4(paridiprubart biosimilar) mAb (BME100611) on HEK293 cells transfected with Human TLR4 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

