

PRODUCT INFORMATION

014763 **Uniprot ID**

GEN1029;GEN-1029;Hx-DR5-01;Hx-DR5-05;iDD 004;iDD-004;iDD004 **Common Name**

Coniugate Unconjugated

DR5;CD262;KILLER;TRICK2;TRICKB;ZTNFR9;TRAILR2;TRICK2A;TRICK2B;TRAIL-R2;KILLER/DR5 Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Formulation & Reconstitution

Host Species IgG type lgG1 Reactivity Human TNFRSF10B Target

Description Anti-TNFRSF10B(benufutamab biosimilar []IgG1) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Background

Usage Research use only

