

PRODUCT INFORMATION

Uniprot ID A0A0B4J1U6

BCD-180;BCD180;BCD 180 **Common Name**

Conjugate Unconjugated

Synonyms TCRBV1S1A1N1;TCRBV9S1

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human **Target** TRBV9

Description Anti-TRBV9(seniprutug biosimilar ☐ IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

