

PRODUCT INFORMATION

Clone ID	OTI2D9 (formerly 2D9)
Target	tRFP
Synonyms	N/A
Host Species	Mouse
Description	Mouse monoclonal turboRFP antibody, clone OTI2D9
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:2000
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated




Figure 1. Anti-tRFP immunoblot analysis of HEK293T over-expression lysates containing recombinant

tRFP, TurboFP602, TurboFP635 fluorescence proteins at 1:2000 dilution.

Address: Wuhan Institute of
Biotechnology B7, Biolake No.666
Gaoxin Road, Wuhan, Hubei, China
Telephone: +1 2409940618(USA)
/+86-18062749453(China)
/+86-400-006-0995(China)

Email: info@dimabio.com
Website: www.dimabio.com



DIMABIO CONFIDENTIAL

