

PRODUCT INFORMATION

P09758 **Uniprot ID**

DS 1062a;DS1062a;DS-1062a;DS-1062;S-1062a;S1062a **Common Name**

Conjugate Unconjugated

EGP1;GP50;M1S1;EGP-1;TROP2;GA7331;GA733-1;TACSTD2 **Synonyms**

Applications N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of Formulation & Reconstitution

reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target** Trop2

Description Anti-Trop2(datopotamab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

temperature.

Research grade biosimilar. Not for use in therapeutic or **Background**

diagnostic procedures for humans or animals.

Usage Research use only

