

## **PRODUCT INFORMATION**

**Uniprot ID** P09758 **Common Name TBC** 

Conjugate Unconjugated

EGP1;GP50;M1S1;EGP-1;TROP2;GA7331;GA733-1;TACSTD2 **Synonyms** 

**Applications** N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions Formulation & Reconstitution

of reconstitution.

**Host Species** N/A lgG1 IgG type Reactivity Human Trop2 **Target** 

**Description** Anti-Trop2(tizetatug biosimilar IgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Research use only Usage

