

PRODUCT INFORMATION

Clone ID	OTI10F11 (formerly 10F11)
Target	tYFP
Synonyms	N/A
Host Species	Mouse
Description	Mouse monoclonal turboYFP antibody, clone OTI10F11
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB IF FC
Recommended Dilutions	WB: 1:2000, IF: 1:50, FC: 1:400
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

