

PRODUCT INFORMATION

Clone ID	OTI7E11 (formerly 7E11)
Target	V5 tag
Synonyms	ANTI-V5; ANTI-V5 TAG; ANTI V5 TAG; ANTI V5-TAG; V5 ANTIBODY; V5 TAG; V5 TAG ANTIBODY
Host Species	Mouse
Description	Mouse monoclonal turboGFP antibody, clone OTI2H8
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:2000
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL






Figure 1. Anti-V5 immunoblot analysis of six different HEK293T over-expression lysates containing different V5 tagged recombinant proteins at 1:2000 dilution.

DIMABIO CONFIDENTIAL

