

PRODUCT INFORMATION

Uniprot ID P08670

ACA 11;ACA-11;ACA11;CLNH 11;CLNH11;CLNH-11 **Common Name**

Conjugate Unconjugated

CTRCT30; epididymis luminal protein **Synonyms**

113;FLJ36605;HEL113

Applications N/A

Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG1 IgG type Reactivity Human VIM **Target**

Description Anti-VIM(pritumumab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

