

PRODUCT INFORMATION

Target	CEACAM5, MUC1
Description	The embedded and sliced tumor sections (4μm section) from Balb/C nu mice implanted with MKN45 cells
Host Cells	MKN45
Applications	IHC
Host mice	Balb/C nu
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage&Shipping	Room temperature
Usage	Research use only




Figure 1. IHC analysis of Balb/C nu MKN45 DiSliceX™ SlideSet section (Cat. No. SLI100006) with Anti-CEACAM5 antibody (DME1012), Rabbit mAb (Cat. No. DME201012) at a concentration of 10 μg/ml. The inset image shows the absence of staining in secondary antibody-only control.

Address: Wuhan Institute of  
Biotechnology B7, Biolake No.666  
Gaoxin Road, Wuhan, Hubei, China  
Telephone: +1 2409940618(USA)  
/+86-18062749453(China)  
/+86-400-006-0995(China)

Email: [info@dimabio.com](mailto:info@dimabio.com)  
Website: [www.dimabio.com](http://www.dimabio.com)



DIMABIO CONFIDENTIAL

