

## **PRODUCT INFORMATION**

**Uniprot ID** P05187

**Common Name** Seagen biosimialr

Conjugate Biotinylated

**Synonyms PLAP** 

**Applications** ELISA, Flow Cyt

Recommended

ELISA 1:5000-10000, Flow Cyt 1:100 **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions.

**Host Species** Humanized

lgG1 IgG type Reactivity Human **Target** ALPP

**Description** Biotinylated Anti-ALPP(SGN-ALPV biosimilar) mAb

**Delivery** 2-3 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

thawing). Lyophilized antibodies are shipped at

ambient temperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

