

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Uniprot ID                   | P05067  |
| Common Name                  | LY-3372993;LY 3372993;LY3372993;N3pG 4;N3pG4;N3PG-IV  |
| Conjugate                    | Biotinylated  |
| Synonyms                     | AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP  |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | IgG1  |
| Reactivity                   | Human   |
| Target                       | APP   |
| Description                  | Biotinylated Anti-APP(remternetug biosimilar) mAb   |
| Delivery                     | 6-8 weeks   |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

