

PRODUCT INFORMATION

Uniprot ID	P21730
Common Name	HIB210,HIB-210;HIB210;MOR044254,MOR-044254;MOR044254;MOR210;MOR-210;MOR 210;TJ210;TJ-210;TJ210;TJ210001;WBP2191;WBP-2191;WBP2191
Conjugate	Biotinylated
Synonyms	C5A;C5AR;C5R1;CD88
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	Human IgG1(K214R,L234A,L235E,G237A,A330S,P331S) – lambda2
Reactivity	Human
Target	C5AR1
Description	Biotinylated Anti-C5AR1(izastobart biosimilar) mAb
Delivery	6-8 weeks
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

