

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Uniprot ID | P20273 |
| Common Name | ACD22-VCMAE, DCDT2980S, FCU2803, RG-7593, RO5541072-000 |
| Conjugate | Biotinylated |
| Synonyms | CD22;SIGLEC2;BL-CAM;SIGLEC-2;Siglec2;SIGLEC2FLJ22814 |
| Applications | ELISA; Flow Cyt |
| Recommended Dilutions | ELISA 1:5000-10000; Flow Cyt 1:100 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Humanized |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD22 |
| Description | Biotinylated Anti-CD22 (pinatuzumab biosimilar) mAb |
| Delivery | 2-3 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

