

PRODUCT INFORMATION

Clone ID 1E9 CD72 **Target Synonyms** Lyb-2 **Host Species** Rabbit

Biotinylated Anti-CD72 antibody(1E9), IgG1 Description

Chimeric mAb

Delivery 2-3 weeks **Uniprot ID** P21854

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human **Applications** Flow Cyt

Recommended

DIMA Disclaimer

Flow Cyt 1/100 **Dilutions**

Purified from cell culture supernatant by affinity **Purification**

chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 Formulation & % – 10% trehalose is added as protectants before Reconstitution lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution. Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Plays a role in B-cell proliferation and **Background**

differentiation.[UniProtKB/Swiss-Prot Function]

Usage Research use only

Biotinylated Conjugate

> All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are

> > Email: info@dimabio.com Website: www.dimabio.com

actively scrutinizing all patent application to

ensure no IP infringement.

