

**PRODUCT INFORMATION**

<b>Common Name</b>	INCA-34176;INCA34176;SNDX-6352;SNDX6352;UCB6352;UCB6352;UCB-6352
<b>Conjugate</b>	Biotinylated
<b>Synonyms</b>	FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R
<b>Applications</b>	N/A
<b>Endotoxin</b>	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG4
<b>Reactivity</b>	Human
<b>Target</b>	CSF1R
<b>Uniprot ID</b>	P07333
<b>Description</b>	Biotinylated Anti-CSF1R(axatilimab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

