Cat. No. BME100489B



## **PRODUCT INFORMATION**

**Uniprot ID** P07333

**Common Name** RG-7155;RG7155;RO-5509554;RO5509554

Conjugate Biotinylated

FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R **Synonyms** 

**Applications** Recommended Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** IgG type lgG1 Reactivity Human CSF1R Target

Biotinylated Anti-CSF1R(emactuzumab biosimilar) mAb Description

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Research use only Usage

