

## **PRODUCT INFORMATION**

**Uniprot ID** P00533

**Common Name** ABX-EGF;E7.6.3;Panitumab;Vectibix

Conjugate

ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2 **Synonyms** 

**Applications** N/A Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

**Host Species** IgG type IgG2 Reactivity Human **Target EGFR** 

Biotinylated Anti-EGFR(panitumumab biosimilar) **Description** 

mAb

6-8 weeks **Delivery** 

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins

Email: info@dimabio.com Website: www.dimabio.com

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background** 

or diagnostic procedures for humans or animals.

Usage Research use only

