

**PRODUCT INFORMATION**

<b>Common Name</b>	IMC001
<b>Conjugate</b>	Biotinylated
<b>Synonyms</b>	GA733-2, M1S2, M4S1, MIC18, TACSTD1, TROP1
<b>Applications</b>	ELISA
<b>Endotoxin</b>	Less than 1.0 EU/μg by the LAL method. For <1 EU/mg requirements, please contact us for customization.
<b>Recommended Dilutions</b>	ELISA 1:5000-10000
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG1 - kappa
<b>Reactivity</b>	Human
<b>Target</b>	EPCAM
<b>Uniprot ID</b>	P16422
<b>Description</b>	Biotinylated Anti-EPCAM(IMC001 biosimilar) mAb
<b>Delivery</b>	2-3 weeks
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

