

PRODUCT INFORMATION

Clone ID	4B2
Target	Exatecan
Synonyms	N/A
Host Species	Rabbit
Description	Biotinylated Anti-Exatecan antibody(4B2), Rabbit mAb
Delivery	In Stock
Uniprot ID	N/A
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	ELISA
Recommended Dilutions	ELISA 1/5000-10000
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Powder
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).
Background	Exatecan (Exa) is a very potent inhibitor of topoisomerase I and anticancer agent. It has been intensively studied as a single agent, a large macromolecular conjugate and as the payload component of antigen-dependent antibody-drug conjugates.
Usage	Research use only
Conjugate	Biotinylated



Biotinylated Anti-Exatecan antibody(4B2), Rabbit mAb ELISA
0.1µg of BSA-Exatecan protein per well




Figure 1. ELISA plate pre-coated by 1 µg/mL (100 µL/well) BSA-Exatecan can Biotinylated Anti-Exatecan antibody(4B2), Rabbit mAb (DME101318B) in a linear range of 16-400 ng/mL.

