

PRODUCT INFORMATION

Uniprot ID	Q99988
Common Name	PF-06946860
Conjugate	Biotinylated
Synonyms	GDF-15;MIC-1;MIC1;NAG-1;PDF;PLAB;PTGFB
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-10000
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	GDF15
Description	Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only



Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb ELISA

0.2 µg of Human GDF15, hFc tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human GDF15 Protein, hFc Tag (PME100628) can bind Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb (BME100090B) in a linear range of 8.0-40 ng/mL.

