

PRODUCT INFORMATION

Uniprot ID P04141

TJ-003234, TJM-2 **Common Name** Conjugate Biotinylated CSF;CSF2 **Synonyms**

Applications ELISA; Flow Cyt

Recommended

ELISA 1:5000-10000; Flow Cyt 1:100 **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species Humanized

lgG1 IgG type Reactivity Human **GM-CSF Target**

Biotinylated Anti-GM-CSF(plonmarlimab biosimilar) mAb **Description**

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

