

**PRODUCT INFORMATION**

<b>Common Name</b>	R-1273;R1273;R 1273;RG-1273;RG 1273;RG1273;RO-4368451;RO4368451;RO 4368451
<b>Conjugate</b>	Biotinylated
<b>Synonyms</b>	NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2)
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	Her2
<b>Uniprot ID</b>	P04626
<b>Description</b>	Biotinylated Anti-Her2(pertuzumab biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage&amp;Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

