

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | YYB101;YYB-101;YYB 101;CLM-101;CLM 101;CLM101;NOV 1105;NOV-1105;NOV1105 |
| Conjugate | Biotinylated |
| Synonyms | SF;HGFB;HPTA;F-TCF;DFNB39 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | HGF |
| Uniprot ID | P14210 |
| Description | Biotinylated Anti-HGF(umikibart biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

