

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Uniprot ID                   | P35225  |
| Common Name                  | DRM 06;DRM-06;DRM06;DRM 6;DRM-6;DRM6;LY 3650150;LY-3650150;LY3650150;MILR1444A;MILR-1444A;MILR 1444;ARG 3637;ARG-3637;ARG3637;RO 5490255;RO-5490255;RO5490255;TNX 650;TNX-650;TXN650  |
| Conjugate                    | Biotinylated  |
| Synonyms                     | P600;IL-13;NC30   |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | IgG4  |
| Reactivity                   | Human   |
| Target                       | IL13  |
| Description                  | Biotinylated Anti-IL13(lebrikizumab biosimilar) mAb   |
| Delivery                     | 6-8 weeks   |
| Storage&Shipping             | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

