

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Uniprot ID | Q16552 |
| Common Name | AIN457 |
| Conjugate | Biotinylated |
| Synonyms | CTLA-8;CTLA8;IL-17;IL-17A;IL17 |
| Applications | ELISA, Flow Cyt |
| Recommended Dilutions | ELISA 1:5000-10000, Flow Cyt 1:100 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Homo sapiens |
| IgG type | Human IgG1 - kappa |
| Reactivity | Human |
| Target | IL17A |
| Description | Biotinylated Anti-IL17A(secukinumab biosimilar) mAb |
| Delivery | 2-3 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

