

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | 15D2, NN-8226, NNC-0109-0012 |
| Conjugate | Biotinylated |
| Synonyms | IL-20;Cytokine Zcyto10 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Homo sapiens |
| IgG type | Human IgG4 - Kappa |
| Reactivity | Human |
| Target | IL20 |
| Uniprot ID | Q9NYY1 |
| Description | Biotinylated Anti-IL20(fletikumab biosimilar) mAb |
| Delivery | 2-3 weeks |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

